singlecellRNA-seq analysis

Pipeline consists of two inputs: metadata file with sample names and condition output from cellranger: barcodes.tsv.gz features.tsv.gz matrix.mtx.gz.

Three R scripts for Integration, dimension reduction and downstream analysis and sh scripts to launch these three files.

One config file

Two python scripts: to create needed output folders

to create nedded folder configuration older with output files from cell ranger. Depending on which kind of the data we received: multiplexed or demultiplexed this Python script will be different.

```
config_seurat.yml dimension_reduction.sh Integration_R.sh
                                                               sample_information.csv
                                                                                                     scRNAseq_do
(singlecell-RNAseq) [maider@login scRNAseq]$ 11
total 60K
drwxr-xr-x 3 maider binf 4.0K Apr 21 12:20
drwxr-xr-x 4 maider binf 88 Apr 20 17:21
drwxr-xr-x 2 maider binf 4.0K Apr 21 14:33 cellranger_run
-rw-r--r- 1 maider binf 1.3K Apr 11 10:14 config_seurat.yml
-rwxr-xr-x 1 maider binf 2.9K Apr 8 15:33 create_proj_str.py
-rw-r--r-- 1 maider binf 139 Apr 8 15:12 dimension_reduction.sh
-rw-r--r-- 1 maider binf 135 Apr 8 15:12 Downstream_Seurat.sh
-rw-r--r-- 1 maider binf 124 Apr 8 15:12 Integration_R.sh
-rwxr-xr-x 1 maider binf 2.0K Apr 21 13:10 link_to_demult_cellranger_output.py
-rw-r--r-- 1 maider binf 42 Apr 14 12:16 sample_information.csv
-rw-r--r-- 1 maider binf 2.0K Apr 20 13:44 scRNASEq_dim_reduction.R
-rw-r--r- 1 maider binf 4.5K Apr 20 13:44 scRNAseq_downstream_analysis.R
-rw-r--r-- 1 maider binf 4.6K Apr 21 16:24 scRNASeq_Integration.R
-rwxr-xr-x 1 maider binf 1.8K Apr 13 14:12 symbolic_cellranger_out.py
(singlecell-RNAseq) [maider@login scRNAseq]$
```

If you are running in the server first:

Ssh maider@ciilogin.c2b2.columbia.edu

And then:

qlogin -l mem=30G

cd /share/data/scRNAseq_test

cp ~/pipelines/scRNAseq/*

The starting working directory will look as follows:

The only input needed to upload at this point from the computer is the csv with the sample metadata.

Sample metadata information will be in csv format as follows. Headers of sample name and condition have to be the same. More columns can be added. Sample name has to be the same as the output folder names in cellranger.

	А	В	С	D	Е	F	
1	orig.ident	Condition					
2	KO1	КО					
3	KO2	ко					
4	WT1	WT					
5	WT2	WT					
6	SERT1	SERT					
7	SERT2	SERT					
8							
9							

scp sample_informatio.csv maider@ciilogin.c2b2.columbia.edu:/path_to_working_directory/

Config.yml has the following variables:

```
data_location: "data"
        cellranger_out: "/share/data/RNA_Seq/10X/Raw_fastq/" #output location for all samples
       cellranger_option: "COUNT" #can be MULTI or COUNT
       metadata_file: "seurat_sample_metadata.csv"
condition_reference: "HC" #condition to use as reference
       quality_location: "Quality_figures/" #output of quality figures
        seurat_object: "Seurat_objects/" #output of Seurat object
10
       seurat_object_name: "Seurat_object.RData" #name of the output Seurat object
       seurat_object_dim: "dim_Seurat_object.RData" #name of the output Seurat obejct with dimensio reduction
11
       integration_folder: "Integration_results/" #outout of downstream integration figures and tables
       DE_folder: "DE_Analysis/" #output of DE analysis
       enrichment_folder: "Enrichment_files/" #output of files to load in GSEA
     #These all are for filtering steps
17
       min_cells: 3
18
       min features: 200
       nFeature_max: 8000
20
       nFeature_min: 50
21
       nCount_max: 60000
       nCount_min: 20
23
       percent_mit: 10
24
       HVF_selection: 2000
25
     #Number of dimensions to use for reduction
27
       npcs: 30
        resolution: 0.3
28
29
       Cell_type_markers: D3D,KLRF1,FCGR3A,CD14,CD1C,IL3RA,IGKC,MYL9,HBB,CD79A
31
      markers_outfile: "Markers.csv" #output name of markers
32
     #Downstream analysis
34
      min.pct: 0.1
        log.fc_threshold_markers: 0.25
35
       log.fc_threshold_DE: 0.25
```

Then we run the python script to create the directories. In the server you will run: python create proj str.py -p /path

```
(singlecell-RNAseq) [maider@login MECFS_test]$ python3 create_proj_str.py -p .
Project Location: .

Checking/Creating sub-directory: ./Quality_figures
=inished creating the./Quality_figures

Checking/Creating sub-directory: ./Seurat_objects
=inished creating the./Seurat_objects

Checking/Creating sub-directory: ./Integration_results
=inished creating the./Integration_results

Checking/Creating sub-directory: ./DE_Analysis
=inished creating the./DE_Analysis

Checking/Creating sub-directory: ./Enrichment_files
=inished creating the./Enrichment_files
=inished creating the./Enrichment_files
```

Then we will copy the output from cellranger. For this we have two python scripts. In the case we received fastq files demultiplexed and we run cellranger count command (go to Commands/Software document) we will use: symbolic_cellranger_out.py.

```
(singlecell-KNAseq) [maider@login sckNAseq]$ (singlecell-RNAseq) [maider@login sckNAseq]$ python symbolic_cellranger_out.py --help usage: This script create the directory structure need to run Seurat R script. Exiting....

Process command line arguments.

optional arguments:
-h, --help show this help message and exit
-cp CELLRANGER_PATH, --cellranger_path CELLRANGER_PATH cellranger output path
-m FILE, --metadata FILE metadata file
```

link_to_demult_cellranger_output.py

(conglectell_RNAseq) [maider@login scRNAseq]\$ python link_to_demult_cellranger_output.py --help
usage: This script create the directory structure need to run Seurat R script. Exiting....

Process command line arguments.

```
optional arguments:
-h, --help show this help message and exit
-cp CELLRANGER_PATH, --cellranger_path CELLRANGER_PATH
cellranger output path
-m FILE, --metadata FILE
metadata file
```

If we receive data multiplexed and we had to run cellranger multi we will use:

In both cases:

```
python symbolic_cellranger_out.py -cp '/share/data/RNA_Seq/10X/Raw_fastq/pbmc_1k_v3_fastqs_2/' -m sample_information.csv
```

```
The data folder format will be as follows:

[(base) mastorkia@maiders-mbp singlecell_MECFS % cd data
[(base) mastorkia@maiders-mbp data % ls

KO1 KO2 SERT1 SERT2 WT1 WT2
[(base) mastorkia@maiders-mbp data % cd KO1
[(base) mastorkia@maiders-mbp KO1 % ls
barcodes.tsv.gz features.tsv.gz matrix.mtx.gz
(base) mastorkia@maiders-mbp KO1 %
```

Once we have all this ready we can launch our scripts.

Go to Seurat_pipeline_documentation